

N-5820P/DRU/LU-101

White LED illumination



This is a LED illumination kit for a simple darkroom unit. It can be attached with a frame (ϕ 28mm pipe) of the unit. Brightness is equivalent to fluorescent lamp 20W.

- LED Color: Main light : daylight color, Auxiliary light: bulb color
- AC Power: 100V/50Hz, 60Hz
- Power consumption: 9W
- Illuminance: 217lx (At the point of 1m from the illumination in a darkroom)
- Total luminous flux: 420lm
- Size: 582×70×44(H)mm
- Weight: 0.96kg
- Power ON/OFF: by a pull code
- Accessories: (1) 5m Extension cord (2) Cable clip ×5pcs

Guide

► Purchase commercial fluorescent lights for replacement.

Specifications

Part Number	N-5820P-D28
-------------	-------------

White LED Light



Specifications

Part Number	LU-101
-------------	--------

Ventilation Fans



Ventilation fans designed to keep out external light. They are two types: an exhaust version (DRU-VF) and a suction version (DRU-IF).

- Size: (W)160×(D)160×(L)193mm (excluding protrusions)
- Weight: 2.2kg
- Power cord: 3m (with switch)
- Power source: AC100V±10% 50/60Hz
- Power consumption: 14W
- Fan suction/exhaust capability: 2.7m³/min. for the motor main unit.

Guide

- Additional processing of black-out curtain is required to install this product in a darkroom. Pipes required for installation come with this product.
- Suction/exhaust capability of 6.2m³/min. type is also available. Contact our Sales Division for more information.

Part Number	Type
DRU-IF	Suction type
DRU-VF	Exhaust type

Head Twist Switch Dust Proof & Water Proof (IP54)

- Pen Clip Included
- Brightness: 18lm
- battery: AAA alkaline batteries × 2 this
- battery Life: About 12 hours
- Accessories: Pen clip, test for battery


LMC/LMT

Levels convenient for installing a vibration isolator, optical laboratory bench, or various other instruments.



- The LMC type is round, and the LMT type is T-shaped.
- The LMC type provides horizontal leveling using the bubble and target. They can be secured on equipment or optical system using their mounting holes.
- The LMT type allows simultaneous leveling in two directions, vertical and horizontal.

Outline Drawing



Specifications

Part Number	External Size [mm]	A [mm]	Sensitivity [°]	Material	Finish	Weight [kg]
LMC-2	φ32	3-φ2.5, PCDφ27	1	Aluminum	Black coating	0.022
LMC-3	φ32	3-φ2.5, PCDφ26	1	Brass	Chrome plated	0.028
LMT-1	45×35	—	1	Aluminum	Black coating	0.027