

Vacuum Pedestal Spacer | VPS-SP

VPS-SP

NEW

RoHS

Catalog Code

W6109




Pedestal spacers for post holder to adjust the height of the optical axis of vacuum lens holders from 12.7mm to 15mm.

- Please insert this spacer between a holder and a post holder (VPST series) and clamp it with M4 screw (VPCB-4-8 or VAGB-4-8, VSC-M4-18EE) .

Outline Drawing

VPS-2.3SP



Guide

- ▶ We deliver the item with specialized vacuum packaging to keep its high clean state from the factory to the destination.
- ▶ We will produce custom products for the requested specifications.

Attention

- ▶ Please open the package in the clean environment. Please do not touch it by your bare hands.

| Specifications | | Primary material : Steel (SUS316L) Finish : None | |
|----------------|----------------------|---|--|
| Part Number | Compatible Holder | Weight | |
| VPS-2.3SP | VMHG-12.7, VLHG-12.7 | 0.01kg | |



Tokyo Head Office
5F, SIGMA KOKI Tokyo Head office 1-19-9,
Midori, Sumida-ku, Tokyo, 130-0021, JAPAN
TEL. +81-3-5638-8228
FAX. +81-3-5638-6550
E-mail international@sigma-koki.com

Osaka Branch
4-9-28 Nishi-Nakajima,
Yodogawa-ku, Osaka 532-0011
TEL. +81-6-6307-4835
FAX. +81-6-6307-4834
E-mail sales.osaka@sigma-koki.com

Kyushu Sales Office
3-17 Hie-machi, Hakata-ku, Fukuoka 812-0014
TEL. +81-92-481-4300
FAX. +81-92-481-4310
E-mail sales.kyushu@sigma-koki.com